



Tria 4.0 Demo Installation Instructions

Thank you for using the Tria Image Processing program.

If you are upgrading from a previous Tria, you may have to remove that installation first. Go to Control Panel/Remove Programs and choose Tria and press uninstall or remove.

Execute the setup.msi program in the installation folder by double clicking on it. This runs the installation and allows you to choose a final destination for the program. You may see a warning about trusting the software. You can trust us.

The most common installation error is a missing Windows driver. This is usually corrected by running Windows Update

Tria uses elements from Microsoft Frameworks 4.8. This is a collection of forms and dialog boxes that simplifies coding, but complicates installation. The installation will notify you if Frameworks is not loaded.

This version is for x64 machines only. Let us know if you need a x86 version.

The TrialImagePack1 contains images used in the processing examples used in the demonstrations in the help file. This image pack can be downloaded from our website.

More Information

Newer Versions and upgrades will be placed on the Quarkt4t website at <http://quarkt4t.com/TriaDemo.html>.

Enjoy!!